



MAN2A1 Antibody, HRP conjugated

Product Code	CSB-PA613697LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16706
Immunogen	Recombinant Human Alpha-mannosidase 2 protein (31-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-mannosidase 2 (EC 3.2.1.114) (Golgi alpha-mannosidase II) (AMan II) (Man II) (Mannosidase alpha class 2A member 1) (Mannosyl-oligosaccharide 1,3-1,6-alpha-mannosidase), MAN2A1, MANA2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MAN2A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.