



DLG1 Antibody, FITC conjugated

Product Code	CSB-PA614264LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12959
Immunogen	Recombinant Human Disks large homolog 1 protein (662-773AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disks large homolog 1 (Synapse-associated protein 97) (SAP-97) (SAP97) (hDlg), DLG1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DLG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.