



TAF1A Antibody, Biotin conjugated

Product Code	CSB-PA614895LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15573
Immunogen	Recombinant Human TATA box-binding protein-associated factor RNA polymerase I subunit A protein (251-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TATA box-binding protein-associated factor RNA polymerase I subunit A (RNA polymerase I-specific TBP-associated factor 48 kDa) (TAFI48) (TATA box-binding protein-associated factor 1A) (TBP-associated factor 1A) (Transcription factor SL1) (Transcription initiation factor SL1/TIF-IB subunit A), TAF1A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TAF1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.