



NEUROD1 Antibody, FITC conjugated

Product Code	CSB-PA615700LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13562
Immunogen	Recombinant Human Neurogenic differentiation factor 1 protein (154-356AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenic differentiation factor 1 (NeuroD) (NeuroD1) (Class A basic helix-loop-helix protein 3) (bHLHa3), NEUROD1, BHLHA3 NEUROD
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NEUROD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.