



NEUROD1 Antibody, Biotin conjugated

Product Code	CSB-PA615700LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13562
Immunogen	Recombinant Human Neurogenic differentiation factor 1 protein (154-356AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenic differentiation factor 1 (NeuroD) (NeuroD1) (Class A basic helix-loop-helix protein 3) (bHLHa3), NEUROD1, BHLHA3 NEUROD
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NEUROD1