




CEBPE Antibody

Product Code	CSB-PA618016LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15744
Immunogen	Recombinant Human CCAAT/enhancer-binding protein epsilon protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CCAAT/enhancer-binding protein epsilon (C/EBP epsilon), CEBPE
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CEBPE

Image



Western Blot

Positive WB detected in: Rat heart tissue

All lanes: CEBPE antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.