




BIRC2 Antibody

Product Code	CSB-PA618777EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13490
Immunogen	Recombinant Human Baculoviral IAP repeat-containing protein 2 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Baculoviral IAP repeat-containing protein 2 (EC 2.3.2.27) (Cellular inhibitor of apoptosis 1) (C-IAP1) (IAP homolog B) (Inhibitor of apoptosis protein 2) (hIAP-2) (hIAP2) (RING finger protein 48) (RING-type E3 ubiquitin transferase BIRC2) (TNFR2-TRAF-signaling complex protein 2), BIRC2, API1 MIHB RNF48
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BIRC2

Image



Western Blot

Positive WB detected in: Rat liver tissue, Rat skeletal muscle tissue, Mouse heart tissue, Mouse lung tissue, Mouse brain tissue

All lanes: BIRC2 antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 70, 65 kDa

Observed band size: 70, 65 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.