



RCAN2 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA618879LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q14206 |
| Immunogen | Recombinant Human Calcipressin-2 protein (88-197AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Calcipressin-2 (Down syndrome candidate region 1-like 1) (Myocyte-enriched calcineurin-interacting protein 2) (MCIP2) (Regulator of calcineurin 2) (Thyroid hormone-responsive protein ZAKI-4), RCAN2, DSCR1L1 ZAKI4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | RCAN2 |