



RCAN2 Antibody, Biotin conjugated

Product Code	CSB-PA618879LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14206
Immunogen	Recombinant Human Calcipressin-2 protein (88-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcipressin-2 (Down syndrome candidate region 1-like 1) (Myocyte-enriched calcineurin-interacting protein 2) (MCIP2) (Regulator of calcineurin 2) (Thyroid hormone-responsive protein ZAKI-4), RCAN2, DSCR1L1 ZAKI4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RCAN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.