



MYBPC2 Antibody, HRP conjugated

Product Code	CSB-PA618884LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14324
Immunogen	Recombinant Human Myosin-binding protein C, fast-type protein (739-1141AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-binding protein C, fast-type (Fast MyBP-C) (C-protein, skeletal muscle fast isoform), MYBPC2, MYBPCF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYBPC2