



MYBPC2 Antibody, FITC conjugated

Product Code	CSB-PA618884LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14324
Immunogen	Recombinant Human Myosin-binding protein C, fast-type protein (739-1141AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-binding protein C, fast-type (Fast MyBP-C) (C-protein, skeletal muscle fast isoform), MYBPC2, MYBPCF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYBPC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.