



MYBPC2 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA618884LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q14324 |
| Immunogen | Recombinant Human Myosin-binding protein C, fast-type protein (739-1141AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Myosin-binding protein C, fast-type (Fast MyBP-C) (C-protein, skeletal muscle fast isoform), MYBPC2, MYBPCF |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | MYBPC2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |