



# RALBP1 Antibody, HRP conjugated

Product Code	CSB-PA618981LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15311
Immunogen	Recombinant Human RalA-binding protein 1 protein (6-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RalA-binding protein 1 (RalBP1) (76 kDa Ral-interacting protein) (Dinitrophenyl S-glutathione ATPase) (DNP-SG ATPase) (Ral-interacting protein 1), RALBP1, RLIP1 RLIP76
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RALBP1