



RALBP1 Antibody, HRP conjugated

Product Code	CSB-PA618981LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15311
Immunogen	Recombinant Human RaIA-binding protein 1 protein (6-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RaIA-binding protein 1 (RaIBP1) (76 kDa RaI-interacting protein) (Dinitrophenyl S-glutathione ATPase) (DNP-SG ATPase) (RaI-interacting protein 1), RALBP1, RLIP1 RLIP76
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RALBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.