



MAPRE2 Antibody, HRP conjugated

Product Code	CSB-PA618987LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15555
Immunogen	Recombinant Human Microtubule-associated protein RP/EB family member 2 protein (159-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated protein RP/EB family member 2 (APC-binding protein EB2) (End-binding protein 2) (EB2), MAPRE2, RP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAPRE2