



MAPRE2 Antibody, Biotin conjugated

Product Code	CSB-PA618987LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15555
Immunogen	Recombinant Human Microtubule-associated protein RP/EB family member 2 protein (159-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated protein RP/EB family member 2 (APC-binding protein EB2) (End-binding protein 2) (EB2), MAPRE2, RP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAPRE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.