



MYLK Antibody, Biotin conjugated

Product Code	CSB-PA618995LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15746
Immunogen	Recombinant Human Myosin light chain kinase, smooth muscle protein (290-398AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin light chain kinase, smooth muscle (MLCK) (smMLCK) (EC 2.7.11.18) (Kinase-related protein) (KRP) (Telokin) [Cleaved into: Myosin light chain kinase, smooth muscle, deglutamylated form], MYLK, MLCK MLCK1 MYLK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYLK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.