



IL10RA Antibody, FITC conjugated

Product Code	CSB-PA621688EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13651
Immunogen	Recombinant Human Interleukin-10 receptor subunit alpha protein (22-235AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-10 receptor subunit alpha (IL-10 receptor subunit alpha) (IL-10R subunit alpha) (IL-10RA) (CDw210a) (Interleukin-10 receptor subunit 1) (IL-10R subunit 1) (IL-10R1) (CD antigen CD210), IL10RA, IL10R
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL10RA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.