



ZMYND11 Antibody, FITC conjugated

Product Code	CSB-PA621870LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15326
Immunogen	Recombinant Human Zinc finger MYND domain-containing protein 11 protein (363-513AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger MYND domain-containing protein 11 (Adenovirus 5 E1A-binding protein) (Bone morphogenetic protein receptor-associated molecule 1) (Protein BS69), ZMYND11, BRAM1 BS69
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZMYND11