



ZMYND11 Antibody, Biotin conjugated

Product Code	CSB-PA621870LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15326
Immunogen	Recombinant Human Zinc finger MYND domain-containing protein 11 protein (363-513AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger MYND domain-containing protein 11 (Adenovirus 5 E1A-binding protein) (Bone morphogenetic protein receptor-associated molecule 1) (Protein BS69), ZMYND11, BRAM1 BS69
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZMYND11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.