



TWF1 Antibody, Biotin conjugated

Product Code	CSB-PA622511LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12792
Immunogen	Recombinant Human Twinfilin-1 protein (199-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Twinfilin-1 (Protein A6) (Protein tyrosine kinase 9), TWF1, PTK9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TWF1