



FOXI1 Antibody, Biotin conjugated

Product Code	CSB-PA622520LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12951
Immunogen	Recombinant Human Forkhead box protein I1 protein (183-343AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein I1 (Forkhead-related protein FKHL10) (Forkhead-related transcription factor 6) (FREAC-6) (Hepatocyte nuclear factor 3 forkhead homolog 3) (HFH-3) (HNF-3/fork-head homolog 3), FOXI1, FKHL10 FREAC6
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FOXI1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.