



TNFAIP1 Antibody, FITC conjugated

Product Code	CSB-PA622686LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13829
Immunogen	Recombinant Human BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 2 protein (233-316AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 2 (hBACURD2) (BTB/POZ domain-containing protein TNFAIP1) (Protein B12) (Tumor necrosis factor, alpha-induced protein 1, endothelial), TNFAIP1, BACURD2 EDP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TNFAIP1