



TNFAIP1 Antibody, Biotin conjugated

Product Code	CSB-PA622686LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13829
Immunogen	Recombinant Human BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 2 protein (233-316AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 2 (hBACURD2) (BTB/POZ domain-containing protein TNFAIP1) (Protein B12) (Tumor necrosis factor, alpha-induced protein 1, endothelial), TNFAIP1, BACURD2 EDP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TNFAIP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.