



IMMT Antibody, Biotin conjugated

Product Code	CSB-PA623101LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16891
Immunogen	Recombinant Human MICOS complex subunit MIC60 protein (76-215AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MICOS complex subunit MIC60 (Cell proliferation-inducing gene 4/52 protein) (Mitochondrial inner membrane protein) (Mitofilin) (p87/89), IMMT, HMP MIC60 MINOS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IMMT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.