



EIF3A Antibody, Biotin conjugated

Product Code	CSB-PA623903LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14152
Immunogen	Recombinant Human Eukaryotic translation initiation factor 3 subunit A protein (1288-1379AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eukaryotic translation initiation factor 3 subunit A (eIF3a) (Eukaryotic translation initiation factor 3 subunit 10) (eIF-3-theta) (eIF3 p167) (eIF3 p180) (eIF3 p185), EIF3A, EIF3S10 KIAA0139
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	EIF3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.