



ALX1 Antibody, FITC conjugated

Product Code	CSB-PA624031LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15699
Immunogen	Recombinant Human ALX homeobox protein 1 protein (10-304AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ALX homeobox protein 1 (Cartilage homeoprotein 1) (CART-1), ALX1, CART1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ALX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.