



dnaK Antibody, HRP conjugated

Product Code	CSB-PA633459HB01EGW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1RGI8
Immunogen	Recombinant Escherichia coli Chaperone protein DnaK protein (1-638AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chaperone protein DnaK (HSP70) (Heat shock 70 kDa protein) (Heat shock protein 70), dnaK
Immunogen Species	Escherichia coli (strain UTI89 / UPEC)
Research Area	Others
Target Names	dnaK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.