



SGO1 Antibody, Biotin conjugated

Product Code	CSB-PA688764LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5FBB7
Immunogen	Recombinant Human Shugoshin 1 protein (260-439AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Shugoshin 1 (Serologically defined breast cancer antigen NY-BR-85) (Shugoshin-like 1), SGO1, SGOL1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SGO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.