



SLC3A1 Antibody, FITC conjugated

Product Code	CSB-PA737352LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07837
Immunogen	Recombinant Human Neutral and basic amino acid transport protein rBAT protein (1-68AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutral and basic amino acid transport protein rBAT (NBAT) (D2h) (Solute carrier family 3 member 1) (b(0,+)-type amino acid transport protein), SLC3A1, RBAT
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC3A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.