



# EIF3M Antibody, HRP conjugated

Product Code	CSB-PA745333LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L2H7
Immunogen	Recombinant Human Eukaryotic translation initiation factor 3 subunit M protein (2-374AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eukaryotic translation initiation factor 3 subunit M (eIF3m) (Fetal lung protein B5) (hFL-B5) (PCI domain-containing protein 1), EIF3M, HFLB5 PCID1
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	EIF3M