




EMB Antibody

Product Code	CSB-PA754239LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PCB8
Immunogen	Recombinant Human Embigin protein (285-327AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Embigin, EMB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EMB

Image



Western Blot

Positive WB detected in: NIH/3T3 whole cell lysate, HepG2 whole cell lysate, Mouse brain tissue, Rat heart tissue


All lanes: EMB antibody at 3μg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 37, 32 kDa

Observed band size: 70 kDa



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA754239LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA754239LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)