




TNPO3 Antibody

Product Code	CSB-PA757558LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5L0
Immunogen	Recombinant Human Transportin-3 protein (382-458AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transportin-3 (Importin-12) (Imp12) (Transportin-SR) (TRN-SR), TNPO3, IPO12
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TNPO3

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA757558LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)