



NPSR1 Antibody, HRP conjugated

Product Code	CSB-PA757818LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6W5P4
Immunogen	Recombinant Human Neuropeptide S receptor protein (1-52AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropeptide S receptor (G-protein coupled receptor 154) (G-protein coupled receptor PGR14) (G-protein coupled receptor for asthma susceptibility), NPSR1, GPR154 GPRA PGR14
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NPSR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.