



# NPSR1 Antibody, Biotin conjugated

Product Code	CSB-PA757818LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6W5P4
Immunogen	Recombinant Human Neuropeptide S receptor protein (1-52AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropeptide S receptor (G-protein coupled receptor 154) (G-protein coupled receptor PGR14) (G-protein coupled receptor for asthma susceptibility), NPSR1, GPR154 GPRA PGR14
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NPSR1