




esxA Antibody

Product Code	CSB-PA763968HA01SKW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GCJ0
Immunogen	Recombinant Staphylococcus aureus ESAT-6 secretion system extracellular protein A protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ESAT-6 secretion system extracellular protein A (Ess extracellular protein A), esxA
Immunogen Species	Staphylococcus aureus
Research Area	Others
Target Names	esxA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: esxA antibody at 2.8μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 27kDa

observed band size: 27 kDa