



esxA Antibody, Biotin conjugated

Product Code	CSB-PA763968HD01SKW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GCJ0
Immunogen	Recombinant Staphylococcus aureus ESAT-6 secretion system extracellular protein A protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ESAT-6 secretion system extracellular protein A (Ess extracellular protein A), esxA
Immunogen Species	Staphylococcus aureus (strain MSSA476)
Research Area	Others
Target Names	esxA