



actbb Antibody

Product Code	CSB-PA768804LA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7ZVF9
Immunogen	Recombinant Zebrafish Actin, cytoplasmic 2 protein (2-375AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin, cytoplasmic 2 (Beta-actin-2) [Cleaved into: Actin, cytoplasmic 2, N-terminally processed], actbb, bactin2
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Signal Transduction
Target Names	actbb
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.