



VASH2 Antibody, HRP conjugated

Product Code	CSB-PA769783LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86V25
Immunogen	Recombinant Human Vasohibin-2 protein (1-355AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vasohibin-2 (Vasohibin-like protein), VASH2, VASHL
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	VASH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.