



FERMT3 Antibody, HRP conjugated

Product Code	CSB-PA771441LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UX7
Immunogen	Recombinant Human Fermitin family homolog 3 protein (388-503AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fermitin family homolog 3 (Kindlin-3) (MIG2-like protein) (Unc-112-related protein 2), FERMT3, KIND3 MIG2B URP2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FERMT3