



FAM19A3 Antibody

Product Code	CSB-PA773794LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z5A8
Immunogen	Recombinant Human Protein FAM19A3 protein (31-169AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM19A3 (Chemokine-like protein TAFA-3), FAM19A3, TAFA3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM19A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.