



TLK2 Antibody, Biotin conjugated

Product Code	CSB-PA801231LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UE8
Immunogen	Recombinant Human Serine/threonine-protein kinase tousled-like 2 protein (23-234AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase tousled-like 2 (EC 2.7.11.1) (HsHPK) (PKU-alpha) (Tousled-like kinase 2), TLK2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TLK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.