



GPIHBP1 Antibody, Biotin conjugated

Product Code	CSB-PA811603LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IV16
Immunogen	Recombinant Human Glycosylphosphatidylinositol-anchored high density lipoprotein-binding protein 1 protein (21-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycosylphosphatidylinositol-anchored high density lipoprotein-binding protein 1 (GPI-HBP1) (GPI-anchored HDL-binding protein 1) (High density lipoprotein-binding protein 1), GPIHBP1, HBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	GPIHBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.