



FOXRED2 Antibody, FITC conjugated

Product Code	CSB-PA816892HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWF2
Immunogen	Recombinant Human FAD-dependent oxidoreductase domain-containing protein 2 protein (501-684AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FAD-dependent oxidoreductase domain-containing protein 2 (Endoplasmic reticulum flavoprotein associated with degradation), FOXRED2, ERFAD
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FOXRED2