



FOXRED2 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA816892HD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8IWF2 |
| Immunogen | Recombinant Human FAD-dependent oxidoreductase domain-containing protein 2 protein (501-684AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | FAD-dependent oxidoreductase domain-containing protein 2 (Endoplasmic reticulum flavoprotein associated with degradation), FOXRED2, ERFAD |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | FOXRED2 |