



FOXRED2 Antibody, Biotin conjugated

Product Code	CSB-PA816892HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWF2
Immunogen	Recombinant Human FAD-dependent oxidoreductase domain-containing protein 2 protein (501-684AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FAD-dependent oxidoreductase domain-containing protein 2 (Endoplasmic reticulum flavoprotein associated with degradation), FOXRED2, ERFAD
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FOXRED2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.