





ARAP3 Antibody

Product Code	CSB-PA819904LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWN8
Immunogen	Recombinant Human Arf-GAP with Rho-GAP domain, ANK repeat and PH domain-containing protein 3 protein (163-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arf-GAP with Rho-GAP domain, ANK repeat and PH domain-containing protein 3 (Centaurin-delta-3) (Cnt-d3), ARAP3, CENTD3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ARAP3

Image



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA819904LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA819904LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.