



ARAP3 Antibody, Biotin conjugated

Product Code	CSB-PA819904LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWN8
Immunogen	Recombinant Human Arf-GAP with Rho-GAP domain, ANK repeat and PH domain-containing protein 3 protein (163-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arf-GAP with Rho-GAP domain, ANK repeat and PH domain-containing protein 3 (Centaurin-delta-3) (Cnt-d3), ARAP3, CENTD3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ARAP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.