




BAS1 Antibody

Product Code	CSB-PA822134HA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96291
Immunogen	Recombinant Arabidopsis thaliana 2-Cys peroxiredoxin BAS1, chloroplastic protein (66-266AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2-Cys peroxiredoxin BAS1, chloroplastic (2-Cys Prx A) (2-Cys peroxiredoxin A) (EC 1.11.1.15) (Thiol-specific antioxidant protein A), BAS1
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	BAS1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: BAS1 antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 28 kDa

Observed band size: 28 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.