



FBXW7 Antibody, Biotin conjugated

Product Code	CSB-PA822163LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969H0
Immunogen	Recombinant Human F-box/WD repeat-containing protein 7 protein (1-707AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/WD repeat-containing protein 7 (Archipelago homolog) (hAgo) (F-box and WD-40 domain-containing protein 7) (F-box protein FBX30) (SEL-10) (hCdc4), FBXW7, FBW7 FBX30 SEL10
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXW7